

# Sprint Assignment 2 Report

By Team D

COEN/ELEC 390  
Department of Electrical & Computer Engineering  
Concordia University

March 11th, 2018

# Table of Contents

List of Tables	2
List of Figures	3
Introduction	3
Product Backlog	3
Sprint 1 Backlog	5
Design Document	6
Android Application Wireframe	6
System Architecture	7
Hardware Architecture	8
Software Architecture	8
Use Cases and Sequence Diagrams	9
Test Document	11
Test Plan 1: Crash Test	11
Test plan 2: Skeleton Reliability	12
Definition of Done	15
Sprint 2 Backlog	16

## List of Tables

<a href="#">Table 1: Product Backlog</a>	5
<a href="#">Table 2: Sprint 1 Backlog</a>	6

## List of Figures

<a href="#">Figure 1: Application Wireframe v2</a>	7
<a href="#">Figure 2: System Architecture v2</a>	8
<a href="#">Figure 3: Reminders Sequence Diagram</a>	10
<a href="#">Figure 4: User Profile Sequence Diagram</a>	11

# Introduction

As it was concluded in Milestone 2, our team decided to continue with the idea of a mobile application for seniors that will remind and warn them to take their medicine on time using notifications

The application will be connected wirelessly to a set of sensors that will be located in a pill box, and will keep notifying the user until the sensors detect activity. The application will also notify a caregiver or a guardian if no activity is detected for a predefined period of time. The application will also educate the patients about the medication taken by providing valuable information to them about the medicines. Those information will contain an illustration of the medication, the recommended intake frequency, the ingredients, the side-effects, what diseases it treats, etc.

We believe that our project will increase the effectiveness of pharmaceutical treatments by resolving the recurring problem of patients that forget to intake medications and by educating them at the same time.

The next step in the development of our mobile application will be the second sprint in which we plan to improve the first version of our software, improve the hardware and establish a connection between both. The software will be the skeleton of our application and will contain basic options such as a login page, a settings page and a reminder scheduler. For the hardware, we plan to create the best possible pillbox that will include two hall-effect sensors. The hardware will be able to detect if the box is open and communicate it to the mobile app through bluetooth.

# Product Backlog

Story ID	Story Title	Card	Story Point	Sprint	Status	Conversation	Confirmation
BK-1	sensor and microcontroller	<i>As a senior person, I want the medication box to be able to detect when I open the pill box, so that I can be reminded to take the medication</i>	8	Sprint 1	Sprint 1	1. The sensors used are hall sensors 2. There will be an Arduino Uno 3. There will be two sensors	1. The sensors can detect a magnetic field 2. The microcontroller can take data from the sensors 3. The sensors are in sync
BK-2	skeleton	<i>As a senior, I want to be able to log into the Medication Taking Reminder app, check the information about my medication, allow my nurse to set the medication intake schedule, and input guardian contact information, so that I may reduce any possibility of missing my medication intake.</i>	13	Sprint 1	Sprint 1	1. It is a prototype of our application. 2. it should have several button for those functions.	1. Does it have a login button? 2. Does it have a check schedule button? 3. Does it have a set medication schedule interface? 4. Does it have set guardian information button? 5. Does each button take the user to a different page?
BK-3	open box detection	<i>As a senior person, I want the medication box to be able to detect when I open the pill box, so that I can get a notification when I miss the time to take my medicine.</i>	8	Sprint 2	Sprint 2	1. Hall sensors are used in the detection 2. A magnet will be used for sensing 3. The box is a common pill box 4. The miss rate for detection needs to be less than 5%	1. Does the sensor successfully detect the "open box" action? 2. Is the miss rate less than 5%?
BK-4	bluetooth connection	<i>As a user of the application, I want the medication box to be able to connect to my smart device so that I can have real-time updates on my device</i>	3	Sprint 2	Sprint 2	1. There will be a Bluetooth slave chip on the sensor module 2. The version is Bluetooth 4.0 3. Sensor data sent via Bluetooth	1. The Bluetooth chip connects to the phone 2. The sensor data is successively transmitted to the phone
BK-5	schedule setting	<i>As a nurse, I want the application to allow me to set a reminder schedule on my patient's cellphone for medication intake, so that I don't need to remind them personally every time for pills.</i>	5	Sprint 2	Sprint 2	1. There will be a schedule planner 2. The guardian will be in charge of schedule planning 3. There will be a form of backup to prevent schedule to be accidentally erased	1. The schedule can be set within 5 minutes 2. Schedule setting is not ambiguous 2. The schedule can be reset 4. The schedule can be updated 5. The schedule has a log to keep track of changes
BK-6	notification setting	<i>As a senior person, I want the application to send me an alert when I don't open the medication box on time, so that I will not miss the chance to take medicine in proper time.</i>	5	Sprint 2	Sprint 2	1. The notification is a message 2. The notification is a sound 3. The notification can be set to have different sounds 4. The user can set the notification frequency and time	1. The notification message needs to include the time to intake, the time of intake, the medication to take and the quantity 2. It is possible to change the notification frequency and times
BK-7	audible notification	<i>As a senior person, I want the notification to be loud enough so that I can be sure not to miss it</i>	5	Sprint 3	Sprint 3	Having this feature will eliminate some of the risks that the notification will not be heard	1. Is the users information only obtainable by specific individuals? (i.e Nurse, Doctor, Guardian, User)
BK-8	notification for guardian	<i>As a guardian, I want the application be able to inform me if my senior does not take their medicine a period time pass the schedule, so that I can contact them and get more detail of why they miss the right time to take pills.</i>	13	Sprint 3	Sprint 3	1. What is the period be we send the notification to the guardian, if the senior does not take their medicine? 2. How can we store the information of guardian? 3. How can the application automatically send the notification to the guardian phone when the senior does not take the pill pass a period of time?	1. Does the notification send when the sensor does not detect a action passing a certain of time? 2. Does the notification successfully send to a particular phone?
BK-9	educative interface	<i>As a senior person, I want to be able to have access to information concerning the medication I take so that I know the benefits, dosage and secondary effects of the medication I'm taking.</i>	20	Future	Not in sprint	Having this feature will allow the user to be better informed on the benefits, dosage and secondary effects of the medication, thus allowing him/her to make better choices and to encourage him/her to take the medication	1. Can the user understand clearly the information provided? 2. Is the information up-to-date with current research and development? 3. Is the information stated in a way which invites the user to read and does not bore them?
BK-10	privacy protection	<i>As a senior, I want to be able to have protected access to the application so that I can preserve my confidentiality</i>	13	Future	Not in sprint	1. What technology are we using to complete this function? 2. Who can access the data?	1. Is the users information only obtainable by specific individuals? (i.e Nurse, Doctor, Guardian, User)
BK-11	multi-box interface	<i>As a senior, I want to be able to monitor my activity when I have more than one medication box so that I do no mix up and can be more organized</i>	20	Future	Not in sprint	1. There will be a page to add, drop and modify box and box content 2. There will be buttons to access these pages on the main page	1. Check if add/drop/modify takes less than 5 minutes to do 2. Check if the buttons get to the different section within 3 clicks
BK-12	multi-sensor detection	<i>As a senior, I want to have a redundant system using multiple sensors to detect my medication activities so that I can be sure I took the medication</i>	40	Future	Not in sprint	1. The suggested sensors are motion sensors, contact sensors, vibration sensors 2. The sensors will complement each other	1. Check sensor data output 2. Check if the controller can handle all the data 3. Check how sensors can complement one another
BK-13	user database	<i>As a guardian, I want to have access to the information concerning the senior from my own mobile device so that I can monitor their activity from any location</i>	20	Future	Not in sprint	1. There will be a database that can be updated through servers 2. The database will hold using information 3. Authorized user can access assigned database	1. Check if updating database is simple 2. Check if authorized users can access theirs but not others 3. Check if servers can process data quickly

**Table 1: Product Backlog**

# Sprint 1 Backlog

Sprint 1 Goal(s):		1. Build sensors and microcontroller circuit and collect data				
		2. Build the skeleton of the application				
Story ID	Task ID	Task Title	Task Description	Ideal Hours	Status	Comments
BK-1	1,1	Brainstorming				
	1.1.1	Brainstorming I	Research about possible sensors that could be possibly used in this application	3 hours	Completed	Motion sensor, vibration sensor, contact sensor, potentiometer, hall effect sensor
	1.1.2	Brainstorming II	Discuss the complexity of implementation associated with the sensor	3 hours	Completed	We finally decided to use the hall effect sensor
	1.1.3	Brainstorming III	Discuss the connectivity to be used between the devices	2 hours	Completed	Bluetooth 4.0 connection
	1,2	Sensor Position	Find out the best position for the sensor	3 hours	Completed	3 at most on the front of the box
	1,3	Obtaining sensors	Research information about different kind of sensors, decide which one is needed, find where the sensors are sold and buy them	2 hours	Completed	Bought them from ABRA electronics, 2 for the moment
	1,4	Circuit design	Design the circuit that will connect the sensors to the microcontroller and allow interfacing with them	4 hours	Completed	A circuit was found on a hobbyist website and it was modified to fit our needs
	1,5	Arduino code	Write the Arduino code that will allow retrieval of data from the sensors	3 hours	Completed	A simple code was written to turn on an led and display closed when the Hall sensor senses a magnet
	1,6	Build/interface sensors	Build the circuit and use the Arduino code to collect data from the sensors	4 hours	Completed	The circuit was built on a bread board and was successfully interfaced with
BK-2	2,1	Login Page	Creation of a login page and links to its related activities	5 hours	Completed	
	2,2	User Information input Activity	user can input their name, guardian's name and guardian's phone number.	5 hours	Completed	
	2,3	User profile Activity	This activity contains four buttons, which are create reminder, view reminder, edit profile, and backward button	5 hours	Completed	
	2,3	Schedule Activity				
	2.3.1	schedule checking activity	Creation of the reading mode of this activity	5 hours	Completed	only allow to read the schedule
	2.3.2	schedule setting activity	Creation of the edit mode of this activity	5 hours	Completed	enables the edition of information
	2,4	User information edit Activity	An activity allows user to change their information.	5 hours	Completed	
	2,5	Action bar	An action bar for those two Schedule activity and edit Activity	5 hours	Completed	
	2,6	Testing				
	2.6.1	Buttons	Test and check that each button in the application meet their respective requirements	5 hours	Completed	
	2.6.2	Schedule Activity	Test and check the schedule activity working well	5 hours	Completed	
	2.6.3	User information edit Activity	Test and check the component in this activity working well	5 hours	Completed	
	2,7	Integration testing	Test the entire application	5 hours	Completed	Involves unit and component testings
	2,8	Defect correction	This is reserved for solving any error that may arise due to task/feature malfunction	5 hours	Completed	
Documentation				6 hours	Completed	
Product Backlog groom				5 hours	Completed	

**Table 2: Sprint 1 Backlog**



# Design Document

## Android Application Wireframe

The wireframe shown in figure 1 is the skeleton of the application. It best interests the needs of seniors since they will occupy the majority of the user base. The user interface will have many big buttons and easy to read text. The creation interfaces only hold 3-4 types of data, that way the user is not overwhelmed with unnecessary options. Also, the user will be able to navigate to any page within two button presses (after the user select screen), this way the user will be able to navigate back to the splash screen with ease. Expect the final application to have similar features to those shown in figure 1.




Figure 1: Application Wireframe v2

## System Architecture

This system architecture will be as shown in figure 2. The hall effect sensors will receive an input from the user and send a signal to the microcontroller which will then send a signal via bluetooth to the device. Once the signal is received, the device will process the signal and shut off the alarm if the correct signal is received. Optionally, the device can send a text message to a guardian if the notification has yet to be disabled for X amount of time. User profiles will be stored in an SQLite table, as well as the reminders which will be held in their own SQLite table per user.




Figure 2: System Architecture v2

## Hardware Architecture

The system will be using two hall effect sensors on either end of the front of the lid of the box. There will be two magnets on the front of the base of the box. The sensors and magnets will line up when the box is closed. Hall effect sensors are used for proximity switching. When the box is closed, the sensor will be next to the magnet, the sensor will respond to the magnetic field and output a certain voltage. Opening the box moves the sensors away from the magnets, thus changing the output voltage. The two different output voltages will be classified as the box being either opened or closed. When a reminder initiates on the device, it will search for an open signal from the microcontroller, when the lid is opened and the signal is received, the application will cease to notify the user of the reminder.

## Software Architecture


The software contains 7 activities. The main activity is the user select, this will be the screen the user sees when they initially launch the app. From the user select screen, the user may proceed into the main functionality of the application by selecting an existing user or create a new user by selecting the create user button. The create user activity currently has three credentials to fill out, username, guardian name, and guardian phone number. When the create/done button is pressed, the credentials will be stored in the user SQLite table, as well as create a new SQLite table for the reminders. When a user is selected, the splash page will appear with 4 options, create reminder, view reminders, edit profile and switch user. The create reminder activity currently has three credentials as well, medication name, day of the week, and time (with future support for an image of the medication). Pressing the create/done button stores the reminder in the current user's reminder table and goes back to the splash page. The view reminders activity will show all existing reminders of the current user in a ListView. Pressing an existing reminder will redirect to the edit reminder activity that will function very similarly to the create reminder activity but will instead update the credentials over creating a new row. There will also be an option to delete a reminder from the ListView in the future, which will remove the reminder's row from the table. The edit user activity functions very similarly to the create user activity but will instead update the credentials over creating a new row. Finally the switch profile button will redirect back to the main activity to select a new or existing user. If a user profile is deleted, the user's table as well as the accompanying reminder table will be deleted.



## Use Cases and Sequence Diagrams

### ***Use Case 1: Creating/Editing & Viewing Reminders (figure 3)***


Starting from the main activity, the user select activity, tapping on user\_btn (User 1, 2, 3...) sends the user to the splash page activity. Tapping on create\_reminder\_btn (Create Reminder) sends the user to the create reminder activity. In this activity there are three credentials to input, a String for the medication name, a String for the day of the week, and an int for the time. There will be future support for a BufferedImage for an image of the medication. When add\_reminder\_btn (Create) is tapped, the information is stored in the SQLite table with INSERT INTO reminder\_table. Alternatively, from the splash page activity, tapping view\_reminder\_btn (View Reminders) sends the user to the view reminders activity. In this activity, the contents of reminder\_table are displayed in ListView format. Tapping a one of the TextViews will prompt the user to edit or delete the reminder. Deleting the reminder calls DELETE FROM reminder\_table and removes the specified row from the database. Editing the reminder sends the user to the same page as it would for creating the reminder but instead of inserting into the table, the table is updated with UPDATE reminder\_table.



***Figure 3: Reminders Sequence Diagram***

### ***Use Case 2: Creating/Editing User Profile (figure 4)***

Starting from the main activity, the user select activity, tapping on create\_user\_btn (Create New User) sends the user to the create user activity. In this activity there are three credentials to input, a String for the users name, a String for the guardian name, and a String for the guardian phone number. When add\_user\_btn (Create) is tapped, the information is stored in the SQLite table with INSERT INTO user\_table. Alternatively, from the user select activity, the user can select a profile to move the the splash page activity. When edit\_user\_btn is tapped, the user is sent to the edit user activity which is the same as the create user activity but instead of inserting into the table, the table is updated with UPDATE user\_table. When delete\_user\_btn is tapped, the user is sent to the delete user activity which is the same as the create user activity but instead of inserting into the table, the table is updated with DELETE FROM user\_table.



***Figure 4: User Profile Sequence Diagram***


# Test Document

## Test Plan 1: Crash Test

*Summary: BK-2, BK-5*

We are trying to see whether the app will crash while switching activities in quick succession for a moderate period of time. If we are constantly switching activities, the application could potentially crash from too many actions happening in a short period of time.

We are trying to see whether the app will crash if we keep adding items to the reminder table. If we too add many reminders to the database, the application may not be able to hold all of them and throw an error or crash.

<b>Test Case 1.1: Stress Test</b>		
<b>Pre-Condition:</b> BK-2: We need a working skeleton to be able to switch between activities.		
Steps	Expected Results	Actual Result
1. Switch to a different activity.	Proper switch through activities	Proper switch through activities
2. Speed up rate of switching to attempt to crash the application.	Application crashes due to too many actions.	Application did not crash.
<b>Pass/Fail:</b>		
		


<b>Test Case 1.2: Overflow Test</b>		
<b>Pre-Condition:</b> BK-2, BK-5: We need a working skeleton and working reminder creation to store reminders		
Steps	Expected Results	Actual Result
1. Navigate to create reminder activity.	Proper switch to correct activity	Proper switch to correct activity
2. Add a reminder to database.	Reminder is added to database.	N/A
3. Continue adding many reminders to the database to crash or generate an error.	Error is generated or crash due to too many insertions to the database.	N/A

**Result:** Not enough data is available to generate a result at the moment.

## Test plan 2: Skeleton Reliability


*Summary: BK-2*


We are going to test the performance of the skeleton. Firstly, we are going to test the login page. The login page should have multiple buttons and they will lead us to different user's profile or user creating Activity. Secondly, we are going to test the user information creation Activity. This Activity has three EditText and one button, and it allows us to input information and create a profile. Thirdly, test the user's profile Activity. This Activity contains four button which are "Create Reminder", "View Existing Reminder", "Edit Profile", and "Switch Profile". Also, each button link to an Activity. Fourthly, Reminder Creation Activity allows use to input the information and set the time for taking medication. Fifthly, Schedule check Activity allows us to check the schedule. Finally, the Edit Profile Activity should allow us to reset the user information.

<b>Test Case 2.1:</b> Login page		
<b>Pre-Condition:</b> Have the application running on a device.		
Steps	Expected Results	Actual Result
1. Click the "Create User" button	Go to User information Activity	Go to User information Activity
2. Click backward icon	Back to login page	Back to login page
3. Click "User 1" button	Go to user' profile Activity	Go to user' profile Activity
<b>Pass/Fail:</b>		
		


<b>Test Case 2.2:</b> User's Information creation Activity		
<b>Pre-Condition:</b> Select the "Create User" button on the login page		
Steps	Expected Results	Actual Result
1. Type in the name	Keyboard show up, and allows input characters	Keyboard show up, and allows input characters
2. Type in the phone number	Number keyboard show up and allows input #	Number keyboard show up and allows input #
3. Click "Create" button	A message shows and return to the login page	A message shows and return to the login page




<b>Pass/Fail:</b>


<b>Test Case 2.3:</b> User's Profile Activity		
<b>Pre-Condition:</b> Select the "User #" button on the login page		
<b>Steps</b>	<b>Expected Results</b>	<b>Actual Result</b>
1. Click "Create Reminder" button	Bring us to Reminder Creation Activity	Bring us to Reminder Creation Activity
2. Click backward icon on the action bar	Back to User's Profile Activity	Back to User's Profile Activity
3. Click "View Existing Reminders" button	Go to Schedule Check Activity	Go to Schedule Check Activity
4. Click backward icon on the action bar	Back to User's Profile Activity	Back to User's Profile Activity
5. Click "Edit Profile" button	Go to User's Information creation Activity	Go to User's Information creation Activity
6. Click backward icon on the action bar	Back to User's Profile Activity	Back to User's Profile Activity
7. Click "Switch Profile" button	Return to login page	Return to login page
<b>Pass/Fail:</b>		
		

<b>Test Case 2.4:</b> Reminder Creation Activity		
<b>Pre-Condition:</b> Select the "Create User" button on the login page, then click "Create Reminder"		
<b>Steps</b>	<b>Expected Results</b>	<b>Actual Result</b>
1. Type in the medication name	Keyboard show up, and allows input characters	Keyboard show up, and allows input characters
2. Select the day from menu	A menu which we can choose between Monday and Sunday	A menu which we can choose between Monday and Sunday
3. Type in the time	Number keyboard show up and	Number keyboard show up and

	allows input #	allows input #
4. Click "Create button"	A message shows and return to User's Profile Activity	A message shows and return to User's Profile Activity
<b>Pass/Fail:</b>		
		

<b>Test Case 2.5:</b> Schedule Check Activity		
<b>Pre-Condition:</b> Select the "Create User" button on the login page, then click "View Existing Reminders"		
<b>Steps</b>	<b>Expected Results</b>	<b>Actual Result</b>
1. Observe this Activity	An empty schedule should be exist.	Blank Activity
2. Click backward icon on the action bar	Back to User's Profile Activity	Back to User's Profile Activit
<b>Pass/Fail:</b>		
		

# Definition of Done

## Definition of Done criteria

After discussing the important criterias that are related to our project, here is the checklist that we decided to use for the validation of each PBI :

- ✓ Produced artifact (code/document) for the PBI
- ✓ Requirements mentioned in PBI are satisfied
- ✓ Acceptance criteria(s) for PBI are satisfied
- ✓ Unit tests performed and passed
- ✓ Integration tests performed and passed
- ✓ Project builds without errors
- ✓ Design Document updated
- ✓ Test Document updated

## Definition of Done validation

Story ID: BK-1	
<b>User Story: Sensor and Microcontroller</b> <i>As a senior person, I want the medication box to be able to detect when I open the pill box, so that I can be reminded to take the medication</i>	
DoD checklist for this PBI:	Status
	Completed
✓ Produced artifact (code/document) for the PBI	Done
✓ Requirements mentioned in PBI are satisfied	Done
✓ Acceptance criteria(s) for PBI are satisfied	Done
✓ Unit tests performed and passed	Done
✓ Integration tests performed and passed	Done
✓ Project builds without errors	Done
✓ Design Document updated	Done
✓ Test Document updated	Done



Story ID: BK-2	
<b>User Story: Skeleton</b> <i>As a senior, i want to be able to log into the Medication Taking Reminder app, check the information about my medication, allow my nurse to set the medication intake schedule, and input guardian contact information, so that I may reduce any possibility of missing my medication intake.</i>	
DoD checklist for this PBI:	Status
	Completed
✓ Produced artifact (code/document) for the PBI	Done
✓ Requirements mentioned in PBI are satisfied	Done
✓ Acceptance criteria(s) for PBI are satisfied	Done
✓ Unit tests performed and passed	Done
✓ Integration tests performed and passed	Done
✓ Project builds without errors	Done
✓ Design Document updated	Done
✓ Test Document updated	Done

## Sprint 2 Backlog

Sprint 2 Goal(s): 1. Successfully detect the "open" action and missing rate should less than 5%. 2. Establish a bluetooth connection between the hardware and the application 3. A notification is sent at a specific time 4. The medicine schedule allows to set, save, and check						
Story ID	Task ID	Task Title	Task Description	Ideal Hours	Status	Comments
BK-3	3,1	Obtaining the box	Find or purchase a suitable box, ideally a commercially available pill box, as our prototype	1 hour	Planned	
	3,2	Box opening action				
	3.2.1	Open definition	Define what is considered as opened and what is considered as closed, this subtask is done by brainstorming	2 hours	Planned	
	3.2.2	Sensor detection algorithm	Make the application able to use information from the sensors	5 hours	Planned	
	3,3	Sensor integration	Create a prototype for the box that will be used for the application by adding the sensors to it	5 hours	Planned	
	3,4	Data collection	Collect data from the sensors for opening and closing the box	2 hours	Planned	
	3,5	Missing rate	Determine the miss rate from the data collected	2 hours	Planned	
	3,6	Reducing the missing rate	If the missing rate is greater than 5%, adjust the sensor to a better position, and calculate the missing rate	5 hours	Planned	
BK-4	4,1	Pill box and smartphone connection				
	4.1.1	Bluetooth integration	Integrate a Bluetooth Chip with the microcontroller system	3 hours	Planned	
	4.1.2	Wireless sensor system	Use a battery (9V) to power the sensor/microcontroller system to make the system independent of the computer's usb port	3 hours	Planned	
	4.1.3	Bluetooth connection	Write the code that will enable the application to receive Bluetooth signals	4 hours	Planned	
	4.1.4	System testing	Test the connection at different distances and in different rooms to verify the connectivity and obtain data	5 hours	Planned	
	4,2	Sensor interfacing	Write the code that will enable the application to use the sensor data received from the sensor module	5 hours	Planned	
BK-5	5,1	Schedule table	use the liveview (SQLite) to create a schedule table to "view existing reminders"	5 hours	Planned	
	5,2	Menu to create reminder	Create a page that allows the user to create, remove and modify reminders	4 hours	Planned	
	5.2.1	Add button	Create a button that will allow creation of a reminder	3 hours	Planned	
	5.2.2	Remove button	Create a button that will allow removal of a reminder	3 hours	Planned	
	5.2.3	Modify button	Create a button that will allow modification of a reminder	3 hours	Planned	



BK-6	6,1	<b>Message notification</b>	Create a notification message that contains time of intake, time to intake, medicine to take and quantity	2 hours	Planned	
	6,2	<b>Notification sending</b>				
	6.2.1	real time function	A function which will get and store the actual time	5 hours	Planned	
	6.2.2	get schedule time	A function which will get and store the time from the schedule	5 hours	Planned	
	6.2.3	compare function 1	A function for comparing the real time and the schedule time	3 hours	Planned	
	6.2.4	check "open" action	A boolean function to check the pill is opened or not	2 hours	Planned	
	6.2.5	<b>notification send</b>	If the pill box does not open, the notification will send. Otherwise, the notification will not send.	3 hours	Planned	
	6.2.6	<b>notification resend</b>	if the pill box does not open 15 mins after the first notification was sent, send the notification again	5 hours	Planned	
Testing and fix defect				10 hours	Planned	
Potential customer feedback				1 hour	Planned	
<b>Documentation (the pdf file)</b>				6 hours	Planned	
Product backlog grooming				5 hours	In Progress	